

FASTCAM NI DAQ
software option

 **User's Manual**

Revision 1.10E

Photron

- The copyright of this manual is held by Photron Limited.
- Product specifications and manual contents can change without advanced notification.
- This manual was created taking every possible measure to ensure the accuracy of its contents. However, if you find a section which is unclear, a mistake, or an omission, please contact Photron Limited using the contact information provided at the end of the manual.
- Photron Limited bears no responsibility for the results of using the product or from following the instructions in this manual.

- The official name of Windows is the Microsoft® Windows® Operating System.
- Microsoft®, Windows®, Windows Vista®, Windows® 7 are registered trademarks or trademarks in the United States and other countries of Microsoft Corporation of the United States.
- Pentium® is a registered trademark of Intel Corporation of the United States.
- National Instruments®, NI®, ni.com® are the registered trademark of National Instruments in the United States (US) and in other countries.
- Other company names and product names listed in this manual are trademarks of their respective companies.
- Product specifications and features can change for the purpose of improvement without notification.

Introduction

Thank you for purchasing your FASTCAM NI DAQ software option.

This manual contains the operating instructions and warnings necessary for using the FASTCAM NI DAQ software option.

Please read the entire manual before using the software.

If any part of this manual is unclear, contact Photron using the contact information printed at the back of this manual.

Manual Notation

The following icons and symbols are used in the explanations in this manual.






Icon/Symbol	Description
 Supplement	This symbol indicates supplementary items to be aware of when using the software.
 Reference	This symbol indicates the location of a reference.
 Important	This symbol indicates content that should always be read.
 Caution	This symbol indicates instructions that should always be followed when using the software, or things to be careful of when using the software.
 MEMO	This symbol indicates a space for you to use for making notes.
“ ”	This symbol is used to indicate the names of items on a screen, references, dialog names, and keyboard keys.
[]	This symbol is used to indicate screen names, button names, and menu names.
< >	This symbol is used to explain operating procedures in diagrams and supplementary items.

Table of Contents

Chapter.1. Overview of the FASTCAM NI DAQ software option	1
1.1. Overview of the FASTCAM NI DAQ software option	2
1.2. Operating environment of the FASTCAM NI DAQ software option	3
1.3. Included accessories of the FASTCAM NI DAQ software option	4
Chapter.2. Installing the FASTCAM NI DAQ software option	5
2.1. Installing the FASTCAM NI DAQ software option	6
2.2. Installing USB driver software for NI DAQ	9
Chapter.3. Using the FASTCAM NI DAQ software option	15
3.1. Recording work flow	16
3.2. Connecting PC and NI DAQ	17
3.3. Connecting FASTCAM and NI DAQ (Manual Trigger)	20
3.4. Connecting FASTCAM and NI DAQ (Level Trigger)	21
3.5. Recording Image	22
3.5.1. Setting the Frame Rate	22
3.5.2. Setting the Shutter Speed	23
3.5.3. Setting the Trigger Mode	24
3.5.4. Trigger Mode Types	27
3.6. Setting the Wave input board bar	30
3.7. Inputting a Rec Start trigger	36
3.8. Graphic Display of Waveform Input Data	37
3.8.1. Settings for Graphic Display of Waveform Data	40
3.8.2. Changing the Channel Name	41
3.8.3. Setting Style of Displayed Graph	42
3.8.4. Setting Display Unit, Gradient and Segment	43
3.8.5. Setting Range of Graphic Display	44
3.8.6. Graph View Window	46
3.8.7. Storing in CSV format	48
Chapter.4. Alert Display	51
Chapter.5. Specifications	53
Chapter.6. Contacting Photron	55

Chapter.1. Overview of the FASTCAM NI DAQ software option

1.1. Overview of the FASTCAM NI DAQ software option

With the FASTCAM NI DAQ software option, it is possible to use “NI USB-6251 BNC” waveform measuring instruments by National Instruments (hereafter referred to as NI DAQ) to simultaneously record analog signals along with images from FASTCAM high speed digital cameras by Photron (hereafter referred to as FASTCAM) and playback graphs and images in synchronization.

The camera may be operated by using a threshold input level to the NI DAQ system, or by START/STOP buttons.

Features of the FASTCAM NI DAQ software option

1. Integrated adjustment of both high speed cameras and waveform measuring instruments

The FASTCAM NI DAQ software option was designed to be used as a plug-in for Photron FASTCAM Viewer (hereafter referred to as PFV), the camera control software for FASTCAM. From this plug-in software, it is possible to perform integrated adjustment of both FASTCAM and NI DAQ.

The FASTCAM NI DAQ software option imports the image recording settings of FASTCAM, and then automatically sets NI DAQ so that it matches with the image recording time of FASTCAM. Furthermore, by saving FASTCAM image data and NI DAQ waveform data at a discretionary location with the same name, it is possible to view data after recording.

2. Includes cables and connectors for simultaneous operation

The product is packaged with the cables and connectors necessary for the simultaneous operation of FASTCAM and NI DAQ. This means that even new users can easily perform setup.

Supplement

- For information regarding applications involving FASTCAM NI DAQ software option, please send your request to Photron or to your nearby Photron dealer.
- If you have any questions regarding NI DAQ, please contact National Instruments.

Caution

- FASTCAM NI DAQ software option is plug-in software and cannot be run in isolation. PFV must be acquired separately.
- Though the “NI USB-6251 BNC” features operation of up to 1.25M samples/s, its operation in FASTCAM NI DAQ software option is limited to 1M samples/s.

1.2. Operating environment of the FASTCAM NI DAQ software option

FASTCAM NI DAQ software option can be run with the following OS environments, high speed cameras, and waveform measuring instruments.

- OS
 - Windows Vista (32-bit or 64-bit version)
 - Windows 7 (32-bit or 64-bit version)
 - Windows 8 (32-bit or 64-bit version)
 - Windows 8.1 (32-bit or 64-bit version)

- High speed camera
 - FASTCAM SA Series
 - FASTCAM MH/MC Series
 - FASTCAM Mini Series
 - FASTCAM Multi

- Waveform measuring instrument
 - NI USB-6251 BNC
(National Instruments)

- PFV
 - Ver.3.6.0 (32-bit/64-bit version) or later

1.3. Included accessories of the FASTCAM NI DAQ software option

The FASTCAM NI DAQ software option includes the following accessories.

- | | | |
|---|--|-----|
| ① | FASTCAM NI DAQ software option installation CD | One |
| ② | FASTCAM NI DAQ software option User's Manual (this manual) | One |
| ③ | Photron FASTCAM Viewer Installation DVD | One |
| ④ | BNC cable (1.5m) | Two |

Chapter.2. Installing the FASTCAM NI DAQ software option

2.1. Installing the FASTCAM NI DAQ software option

Install the FASTCAM NI DAQ software option using the following procedure.


Caution

- Before installing the FASTCAM NI DAQ software option, it is necessary to install PFV.
- FASTCAM NI DAQ software option supports PFV Ver.3.6.0 or later.

Reference

- For the PFV installation procedure, see "1.2 Setup" of the *PHOTRON FASTCAM Viewer User Manual*.


- ① Place the FASTCAM NI DAQ software option installation CD in the CD (or DVD) drive, and then open the CD (or DVD) drive from the computer.
- ② On the CD (or DVD) drive, open the "Setup32" folder, and then double click the file "Setup.exe". The setup program starts, and then the following dialog box appears. Click [Next].



Supplement


- For 64-bit OS, start "Setup.exe" in the "Setup64" folder.

- ③ The user license agreement appears. Confirm the contents, and then click [Yes].



- ④ Specify the installation location. The default setting for the installation location is as follows:
C: \Program Files\PHOTRON\Photron FASTCAM Viewer 3


If necessary, click [Browse] and specify another location. Click [Next].



! Caution

- For directory for installation, be sure to specify a folder where Photron FASTCAM Viewer 3 has been installed.

- ⑤ Installation starts.
When installation is finished, the following screen appears. Click [Finish].




2.2. Installing USB driver software for NI DAQ

To use FASTCAM NI DAQ software option to control NI DAQ, USB driver software must be installed in advance.

Install the USB driver software using the following procedure.

- ① The NI-DAQmx X.X.X DVD comes with NI DAQ. Place the DVD in the DVD drive.



- ② Setup program starts up and the below dialog appears. Click [Install NI-DAQmx X.X.X]




Supplement


- [X.X.X] in the above description indicates the version number.

- ③ Specify the installation location. The default setting for the installation location is as follows:
C: \Program Files\National Instruments\


If necessary, click [Browse] and specify another location.
Click [Next].




- ④ Select an install option.
Select [Typical] and click [Next].




- ⑤ Confirmation process of product information begins.
If your PC is connected to an outside network, there may be a case where certain relevant information is displayed.




- ⑥ License agreement window is displayed. Check the content and select [I accept the above 4 License Agreement(s)] and click [Next].




- ⑦ License agreement widow is displayed again. Check the content and select [I accept the above 2 Licenses Agreement(s)] and click [next].




- ⑧ Trust confirmation window toward the driver software is displayed. Check the content and select [Always trust software from National Instruments Corporation] and click [Next].




- ⑨ A group of software applications being installed are displayed. Click [Next].




- ⑩ Installation of driver software begins. It may take around 30 minutes to install (depending on the PC's performance).



- ⑪ Below window appears when driver software installation is finished.



- ⑫ Below window appears prompting reboot of the PC. Click [Restart] to restart the PC.



! Caution

- To use FASTCAM NI DAQ software option to control NI DAQ, USB driver software included on the NI-DAQmx X.X.X installation DVD or the latest USB driver software downloaded from the National instruments web site must be used. The latest USB driver software for NI DAQ and installation manuals can be acquired from the following website.

<http://www.ni.com/support/>

Chapter.3. Using the FASTCAM NI DAQ software option

3.1. Recording work flow


Use the following procedure to use the FASTCAM NI DAQ software option to record and save FASTCAM image data and NI DAQ waveform data or playback graphs and images in synchronization.

- ① Connect PC and NI DAQ using the specified connection method.
(See 3.2 Connecting PC and NI DAQ)
- ② Connect FASTCAM and NI DAQ using the specified connection method.
(See 3.3 Connecting FASTCAM and NI DAQ (Manual Trigger),
or 3.4 Connecting FASTCAM and NI DAQ (Level Trigger))
- ③ Using PFV, set the various image recording settings, including FASTCAM frame rate, shutter speed, and trigger mode.
(See 3.5 Recording Image)
- ④ Using the Synchronized Navigator Window of the FASTCAM NI DAQ software option, set NI DAQ to match the FASTCAM image recording settings.
(See 3.6 Setting the Wave input board bar)
- ⑤ Recording by FASTCAM and NI DAQ.
(See 3.7 Inputting a Rec Start trigger)
- ⑥ After completion of recording by FASTCAM and NI DAQ, the image data recorded by FASTCAM can be played back using PFV
(See 3.8 Graphic Display of Waveform Input Data)

3.2. Connecting PC and NI DAQ

Before using the FASTCAM NI DAQ software option, connect PC and NI DAQ using the following method.

- ① Prepare the NI USB-6251 USB main unit. Connecting terminals are located on the top.



- ② Connect the power adapter cable to [11-30VDC, 20W] terminal on the main unit. Turn the outer shell clockwise to lock.



- ③ Plug the male connector (B type) on the attached USB cable into female USB receptacle (B type) on the main unit.




- ④ Plug the male connector (A type) on the attached USB cable into female USB receptacle (A type) on PC.




- ⑤ Connect the AC adapter power cable to NI DAQ unit. Plug the AC cord from the adapter to AC outlet.



- ⑥ Turn on power switch of NI DAQ unit. Green LED lights on upper right corner of the main unit.




- ⑦ The NI DAQ should be detected by PC at this point. Installation of the NI DAQ will automatically be completed if the driver software is properly installed by “2.2 Installing USB driver software for NI DAQ”



3.3. Connecting FASTCAM and NI DAQ (Manual Trigger)

This section discusses settings for manual shooting using the REC button on the PFV. Connect between your FASTCAM and NI DAQ as shown below before using the FASTCAM NI DAQ software option.




Caution

- When connecting FASTCAM and PC to the NI DAQ, use either the cables included with FASTCAM and NI DAQ or the BNC and USB cables, included with the FASTCAM NI DAQ software option.
- For the above Gigabit LAN connection, use network terminals compatible with "1000Base-T (Gigabit Ethernet)" and cables compatible with "Enhanced category 5 (CAT5e)" or higher only.
- START, CENTER, END and MANUAL trigger modes can only be used for manual shooting. No other modes are available for use.
- To use a manual trigger, a pre-trigger period of 2 frames or more is needed.

3.4. Connecting FASTCAM and NI DAQ (Level Trigger)

This section discusses settings for shooting that uses the input level difference of waveform data acquired by NI DAQ, higher or lower than a predetermined voltage level, as a trigger. Connect between your FASTCAM and NI DAQ as shown below before using the FASTCAM NI DAQ software option.




Caution

- When connecting FASTCAM and PC to the NI DAQ, use either the cables included with FASTCAM and NI DAQ, or the BNC and USB cables, included with the FASTCAM NI DAQ software option.
- For the above Gigabit LAN connection, use network terminals compatible with “1000Base-T (Gigabit Ethernet)” and cables compatible with “Enhanced category 5 (CAT5e)” or higher only.
- For a Level Trigger shooting, use a signal input to “APF10” or “Dev1/ai0~15” terminal to determine the level difference. However, if a signal to one of “Dev1/ai0~15” terminals is used, no other terminals cannot be used. So, the use of “APF10” terminal is recommended.
- CENTER, END and MANUAL trigger modes can only be used for Level trigger shooting. No other modes are available for use.
- To use input level triggering, a pre-trigger period of 2 frames or more is needed.


3.5. Recording Image

3.5.1. Setting the Frame Rate

The conditions for the shot, such as the frame rate, shutter, and resolution, must be set before recording the subject. Shot conditions for the camera are set from the “Camera” tab on the control panel.



To set the frame rate or the resolution, display the “Camera” tab on the control panel. The frame rate setting can be selected from the list by clicking the [Frame Rate] button. The resolution setting can be selected from the list by clicking the [Resolution] button. The list displayed varies depending on camera model.




3.5.2. Setting the Shutter Speed

Display the “Camera” tab on the control panel and click on the [Shutter] button, a list of shutter speeds that can be set on the connected camera will be displayed.

Select the desired speed from the list.

The list displayed differs depending on the camera model.




3.5.3. Setting the Trigger Mode

Depending on the trigger mode set, the timing of the frames, when they are recorded and the number of frames recorded on the camera when the trigger signal is input will differ.

Caution

- Trigger modes that FASTCAM NI DAQ software option supports are START, CENTER, END and MANUAL - a total of four trigger modes.
- When using a level triggering method, START trigger mode is NOT available for use.
- When using CENTER, END or MANUAL mode, an error is issued if a pre-trigger period is not used completely.


■ Relationship Diagram for the Trigger Modes and Recorded Frames



- Record Button


Use the [Record] button on the control panel to start recording.

The name of the [Record] button changes according to the camera's state. The button name indicates the current state




■ Endless Record

When the images directly before the trigger signal input are recorded (recording with center mode, manual mode, end mode), endless recording is performed.
 When the memory capacity is surpassed, the frames are overwritten.



3.5.4. Trigger Mode Types



! Caution

- Trigger modes available and displayed are dependent on the camera model.
- An error will result if a trigger mode that is NOT supported is chosen for shooting.


■ Start Mode

In start mode, recording starts at the same time the trigger is input.

This mode is suitable for recording a high-speed phenomenon where the events start point is known in advance.


This mode records until the maximum recordable frame count (camera memory capacity).

Example: For a camera that can record 2 seconds in total
2 seconds of high-speed images are recorded directly after the [Trigger In]
button is clicked.




■ Center Mode

In center mode, the content directly before and after when the trigger is input is recorded. This mode is suitable for recording before and after an important moment. This mode records until the maximum recordable frame count (camera memory capacity). The before and after frame count are equal.



■ End Mode


In end mode, the content directly before the trigger is input is recorded. This mode is suitable for recording an important phenomenon for which it is difficult to predict when the event will begin and end. This mode records until the maximum recordable frame count (camera memory capacity).



- **Manual Mode**


Similar to center mode, this mode records the content directly before and after when the trigger is input.

In center mode, the proportion of the recording before and after the trigger input is equal, but in manual mode, the proportion (seconds or frames) recorded before and after the important moment can be set to the desired value.




When recording in manual mode, the settings below must be made in prior to beginning the recording.

1. Click the [Trigger Mode] button and select "Manual".



2. Specify the frame count to record before and after the trigger input. Drag the "Total Frames" slider or enter the frame count directly. When the setting is finished click the [Close] button.



3.6. Setting the Wave input board bar

Set up the PFV waveform port bar.


Before starting PFV, make sure that FASTCAM, NI DAQ and control PC are properly connected and powered up.

- ① Press [START] button, then [All programs] – [Photron FASTCAM Viewer 3] and [PFV Ver.3].

Or, double click  button on the desktop.

- ② PFV starts up.

Input waveform board bar is displayed if FASTCAM NI DAQ software option is installed properly.



Supplement

- If the waveform input board bar is not displayed, click [View] – [Waveform Input Board Bar] on the menu.

- ③ Click the [Open] button on the waveform input board.




- ④ With “Select Device”, select the “NIDAQ”.
Select a board to use from the list displayed in “Device Name” and click the [OK] button.



 **Supplement**

- The name to be displayed in “Device name” should be the same as what is shown in “Measurement & Automation”. To change the device name, see the relevant section in NI DAQ manual.

- ⑤ “Analog Board Setup” dialog is displayed.
Make necessary settings on NI DAQ.



⑥ Settings for the waveform input board are shown below:

The screenshot shows the 'Analog Board Setup' dialog box with several sections highlighted by red boxes and numbered callouts:

- 1**: Sampling section containing 'Samples / frame' (10), 'Sampling Frequency' (10000.00 [Hz]), 'Clock Type' (Internal), 'Number of Sampling' (27270), and 'Input Method' (Differential).
- 2**: Alert section with an empty text field.
- 3**: Channel section containing a table of channels and an 'Input Range' dropdown.
- 4**: Option section with checkboxes for 'Use RecPos Signal', 'Auto Setup', and 'Use Analog trigger to trigger the camera'.
- 5**: External Out section with checkboxes for 'Trigger Out (PFI5)' and 'Sync Out (PFI4)'.
- 6**: External In Terminal section with 'Stop Trigger In : PFI0'.
- 7**: Trigger Condition section with 'Start' (Software, Level 3) and 'Stop' (External Down, Level 3) settings.

Channel	Channel Name	Live [V]
<input checked="" type="checkbox"/>	Dev 1/ai0	Dev 1/ai0
<input type="checkbox"/>	Dev 1/ai1	Dev 1/ai1
<input type="checkbox"/>	Dev 1/ai2	Dev 1/ai2
<input type="checkbox"/>	Dev 1/ai3	Dev 1/ai3
<input type="checkbox"/>	Dev 1/ai4	Dev 1/ai4
<input type="checkbox"/>	Dev 1/ai5	Dev 1/ai5
<input type="checkbox"/>	Dev 1/ai6	Dev 1/ai6
<input type="checkbox"/>	Dev 1/ai7	Dev 1/ai7

Input Range: -10 to 10V

1. Sampling
Sets options related to sampling.

Item	Description
Samples / frame	Select from the list how many points of data to acquire in 1 frame. The sampling frequency is set according to this value and the camera's frame rate. When setting a desired sampling frequency, select "Off".
Clock Type	Select the clock type. Internal : Internal clock, operates with waveform input board's internal clock. External : External clock, operates with the rise or fall of an externally input clock signal. When using an external clock, you must also set the sampling frequency of the signal to input. When setting the external clock's sampling frequency matched to the camera's V-Sync output frequency, check "Sync Camera V-Sync". For normal operation, it is recommended to select "External" and check "Sync with camera V-sync".
Number of Sampling	Sets the number of samplings. If "Auto Setup" is checked, this setting is automatically calculated from the sampling frequency, frame rate, and recording frame count. If an alert is displayed, adjust manually.
Input Method	Indicates the connection method for analog input signal. Choose one from [Referenced Single-Ended], [Non-referenced Single-Ended] and [Differential].

 **Caution**

- If the waveform input board is set to a sampling frequency that is not possible in the specification, it is automatically set to an approximate value.

 **Supplement**

- If "Samples / frame" is set to "Off", sampling frequencies that fulfill the below conditions are displayed in the "Sampling Frequency" list.
 - Integral multiple of the frame rate
 - Input board specification

2. Alert
An alert is issued here when a problematic setting has been made.
Use the displayed information for problem solving and/or asking Photron for help.
3. Channel
Selects channels to use.
Beware of the decrease in measurable sampling speed as the number of channels increases.
4. Option
Sets a collaborative method with the camera.

Item	Description
Use RecPos Signal	Check this item to use the camera start signal, issued by camera, to control the measurement sequence of NI DAQ system. This item is mutually exclusive against "Use Analog trigger to trigger the camera".
Auto Setup	Most of other settings are automatically made by checking this item. Basically, it is recommended to keep this item always checked.
Use Analog trigger to trigger the camera	Check this item if you wish to have the camera triggered by sending a start signal issued at detection of a trigger based on input voltage to NI DAQ. This item is mutually exclusive against "Use RecPos Signal".

5. External Output
Use this item when you wish to have various signals output by NI DAQ.
If "Trigger camera with analog trigger signal" is chosen at "Option", "Trigger output (PF15)" will be checked automatically.
6. External Input
Indicates which signal should best be input to which terminal.
7. Trigger Conditions
Settings to have camera work in collaboration with NI DAQ.
When both "Auto setup" and "Use RecPos Signal" are checked, no display is made because all items are automatically set.

! Caution

- To use "Use Analog trigger to trigger the camera", be sure to make a test shooting to confirm proper operation. This is vitally important because the camera fails to start shooting if there is even the slightest error in setting.

Item	Description	
Start	Sets the recording start and stop conditions.	
	Software	Input the trigger in conjunction with the PFV Record button
	External Down	External trigger falls
	External Up	External trigger rises
Stop	Sample Num	When the specified number of samples is obtained (stop trigger only)
	Level Up	When the voltage input to the channel passed the voltage level set in [Level] in the rising direction
	Level Down	When the voltage input to the channel passed the voltage level set in [Level] in the falling direction
Channel	<p>When “Level Up” or “Level Down” is specified for trigger condition (“Start” or “Stop”), a channel is assigned to accept input analog signal, which is to be used as analog trigger signal.</p> <p>For a Level Trigger shooting, either “APF10” or “Dev1/ai0~15” may be specified. However, the use of “APF10” terminal is recommended basically.</p>	
Post Trigger Deley	<p>Enter the delay count.</p> <p>From the input of the trigger, start/end the actual sampling of data delayed by the amount of the delay count entered here.</p>	


Caution

- When “Auto Setup” is checked, the settings change in conjunction with changes in the camera settings.
- “APF10” or terminals “Dev1/ai0~15” may be specified as a channel for input level detection. However, if one of “Dev1/ai0~15” terminals is specified, ONLY that specified channel can be used for measuring input voltage.


3.7. Inputting a Rec Start trigger

To start a recording on the FASTCAM camera in collaboration with NI DAQ, there are two different ways.

- ① To start measurement on NI DAQ at a start trigger signal issued by the camera by pressing the “Record” button on the PFV.
The whole operation starts going on when the “Record” button is pressed on PFV. The NI DAQ system automatically starts and ends its measurement cycle.




- ② To trigger camera for shooting at an input level trigger issued by the NI DAQ.
A trigger signal is sent to the camera from the waveform input board appropriately set up to certain conditions at “Trigger setting”.
The camera should be set in advance for [Trigger In] or for [Endless Rec] by REC button.



3.8. Graphic Display of Waveform Input Data

This step sets the camera to graphically display external data at a recording.


- ① After recording, have the [Camera] tab displayed on the control panel and, within the [Camera] tab, select the camera that has recorded the waveform input data to be replayed.




- ② Select and check [Memory] in the [Display] field.




- ③ Click the [Graph] button. Note, however, if waveform input data has not been recorded during the recording, this button is grayed out.




- ③ The desired data is downloaded and [Graph Setup] dialog box is displayed. The camera name or waveform input board name that captured the waveform input data is displayed in "Device". Check Analog to select the recorded data for display.



- ⑤ Select channels to display.
Check the desired channels in the [Select] column and click the [OK] button.




- ⑥ Graph is displayed.
• An example of graphically displayed analog waveform data.



3.8.1. Settings for Graphic Display of Waveform Data

After download of waveform input data, or by right-clicking on the graph and clicking [Graph Setup], the [Graph Setup] dialog box appears.


In the [Graph Setup] dialog box, parameters such as the width, color and gradient for the graph can be set.




Item	Description
Select	Sets Hide/Unhide for the recorded data.
Device	Shows the name of the camera that has recorded waveform data.
Channel	Shows the channel name. The name can be changed by double clicking on it.
Line style	Shows the style of lines to display graph - color, type, width, etc. of lines. A right click can change line styles.
Y-Unit	Shows the unit of displayed data. A right click can change units.
Slope	Shows the slope of graph of displayed data. A right click can change slope values.
Y-Intercept	Shows the Y-intercept of displayed data. A right click can change intercept values.

3.8.2. Changing the Channel Name

- ① Select the channel to change and click the [Channel] button.




- ② Enter the desired name in "Channel Name" and click the [OK] button.
The channel name is reset to the original channel name (displayed in "Channel") by clicking the [Reset] button.
All the channel names are reset by clicking the [Reset All] button.




3.8.3. Setting Style of Displayed Graph

- ① Click on the channel where setting should take place.




- ② Click [Line style] and select a color, style and width from the displayed list. Setting can also be made by right clicking on the [Line style] column of selected channels.




3.8.4. Setting Display Unit, Gradient and Segment

- ① Click the [Conversion factor] button.



- ② Set a unit from the [Y Unit] list and click the [OK] button. Select [user] and set any value.



! Caution

- Slope and Y-Intercept can be set to each channel independently.
A unit cannot be set independently on each channel. One same unit is set and reflected on all channels.


3.8.5. Setting Range of Graphic Display

◆ Axis (scale or calibrations)

A range of display can be set.


Check the [Auto Scaling] checkbox of [T] field, and the time axis range of the displayed graph is interlocked with the playback range on the play control panel.

Check the [Auto Scaling] checkbox of [Y] field, the Y-axis range is automatically controlled and readjusted so that all the data in the Y direction can be displayed.



◆ T-Unit


The unit for the time axis set is set.



◆ Display of current bar

Check the [Center] checkbox and the current bar is fixed in the center of the display screen. In usual playback operations, the image display screen and the current bar on the graph are displayed interlocked with each other.


When “Center” is checked, only the maximum value of the time axis can be changed (the minimum value is fixed at 0).



◆ Legends

When [On/Off] is checked, the legend window is displayed.


The information provided by legend covers the maximum and minimum values and the current value.







3.8.6. Graph View Window

◆ Graph View Window Button Functions

You can set the following with the buttons on the graph view window.




Button	Button Name/Function
	Area Zoom After clicking this button, drag the cursor drawing a rectangle over the portion of the graph area to display and the specified area is enlarged.
	Auto-Scale Displays the entire waveform.
	Frame Rewind Rewinds the frame by the data point unit (index unit) when clicked.
	Frame Advance Advances the frame by the data point unit (index unit) when clicked.

◆ Adjustment display area


A mouse drag in the horizontal direction on the graph can move the position of the current bar.

A mouse drag in the vertical or horizontal direction along the Y-axis or T-axis, respectively, changes the display area.




3.8.7. Storing in CSV format


- ① Graphically display the waveform input data that you wish to store in the PC.
- ② Right click on the graph and select [Export Data].




- ③ The "Export data" dialog box is displayed. From "Camera List", select the camera or waveform input board that captured the data to save.



- ④ Specify a frame range to be stored using the "Output frame range" window.




- ⑤ Select an output format by "Output data type", check data type checkboxes and click the [Export] button.



Reference

- For the CSV output format, refer to "8.9. CSV File Structure" in the PFV User's Manual.

- ⑥ Specify the file name and the storage to store the file in, and click the [Save] button.



- ⑦ After export is done, click the [OK] button.
Click the [Close] button of the "Graph data output" dialog box to close it.

Chapter.4. Alert Display

When there is an error in a waveform input board setting, the [Alert] button on the waveform input board bar is displayed in red.

If you click the [Alert] button, a message box is displayed and you can confirm the content of the alert. This error is also displayed in "Alert" on the "Analog Board Setup" dialog box.

When "Auto Setup" is checked, an alert may be displayed when the frame rate or resolution is changed.



Error Display	Error Content
Frequency Over	The sampling frequency exceeds the range that can be set.
Samples / frame not integer error	The sampling frequency is not an integral multiple of the frame rate.
Frequency change Error	The intended sampling frequency was not set. The sampling frequency has been set with the nearest possible value.
Number of Sampling Over	The number of samples exceeds the value that can be set.
Camera Trigger mode Error	The present trigger mode cannot be synchronized with the waveform input board.
Channel Error	The number of channels that can be set has been exceeded.

Chapter.5. Specifications

Product specification serves as the following.

Item	Specification
Analog Input	
Isolated specification	Unisolated
Type	Single-ended , Referenced Single-ended , Differential
Input channels	8ch
Input range	$\pm 10V$ 、 $\pm 5V$ 、 $\pm 2V$ 、 $\pm 1V$ 、 $\pm 0.5V$ 、 $\pm 0.2V$ 、 $\pm 0.1V$
Input impedance	10G Ω
Resolution	16bit
Sampling rate	1.00 MS/s (It divides at two or more times of ch use)
General Specification	
Size	286 (W) x 170 (H) x 69 (D) mm Excluding protrusions
Weight	42.3oz (1.2kg)
Operating environment	
Operating temperature	0 to 45°C
Storage temperature	-20 to 70°C
Humidity	10 to 90% (No Condensation)
Pollution degree	2
Overvoltage category	I

Chapter.6. Contacting Photron

For inquires related to PFV, contact Photron at the contact information listed below. Additionally, the following items will be verified when inquiring, so please prepare them in advance.

Items Verified	Concrete Example
Contact Information	Company, school or organization name, customer contact name, contact phone number, contact e-mail.
Product Name	The Photron FASTCAM Viewer version number and the high-speed camera name. For the version number, check the version information.
Condition of the system and what is known about it.	

Contact Information	
In Americas and Antipodes	<p>PHOTRON USA, INC. 9520 Padgett Street, Suite 110 San Diego, CA 92126-4426, USA Phone : 800-585-2129 or 858-684-3555 Fax : 858-684-3558 E-mail : image@photron.com www.photron.com</p>
In Europe, Africa and India	<p>PHOTRON EUROPE LIMITED The Barn, Bottom Road, West Wycombe, Buckinghamshire, HP14 4BS, U.K. Phone : +44(0) 1494 48 1011 Fax : +44(0) 1494 48 7011 E-mail : image@photron.com www.photron.com</p>
In other areas	<p>PHOTRON LIMITED 21F, Jimbocho Mitsui Bldg., 1-105 Kanda Jimbocho, Chiyoda-Ku, Tokyo 101-0051 Phone : +81 3 3518 6271 Fax : +81 3 3518 6279 E-mail : image@photron.co.jp www.photron.co.jp</p>

Photron *FASTCAM Viewer*

for High Speed Digital Imaging

FASTCAM NI DAQ software option
User's Manual Revision 1.10E

Publication Date July, 2015

Publisher PHOTRON LIMITED

21F, Jimbocho Mitsui Bldg.,

1-105 Kanda Jimbocho, Chiyoda-Ku, Tokyo 101-0051